

# EBP09420 - IMMUNIZING PEPTIDE

**SKU:** EBP09420



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_003033.3

**Immunogen** Peptide with sequence C-EQPQAGSSTSKE, from the C Terminus of the protein sequence according to NP\_003033.3.

**Product Comments** The peptide was used in the production of EBP09420 - Goat Anti-GAT1 / SLC6A1 Antibody. It is currently untested in blocking applications.

**Peptide Sequence** C-EQPQAGSSTSKE

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png)

## DOCUMENTS

- [Data Sheet](#)