

# EBP09426 - IMMUNIZING PEPTIDE

**SKU:** EBP09426



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	
<b>Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_031392.3; NP_001007268.1; NP_001182570.1
<b>Immunogen</b>	Peptide with sequence C-EKKKDTPKQHGT, from the Internal region of the protein sequence according to NP_031392.3; NP_001007268.1; NP_001182570.1.
<b>Product</b>	The peptide was used in the production of EB09426 - Goat Anti-PLA2R1 Antibody. It is currently untested in
<b>Comments</b>	blocking applications.
<b>Peptide Sequence</b>	C-EKKKDTPKQHGT
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)