

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09493 - IMMUNIZING PEPTIDE**

**SKU:** EBP09493

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_002200.2; NP\_001107852.1 ID

 $\label{eq:local_$ 

The peptide was used in the production of EB09493 - Goat Anti-CD11A / integrin alpha L Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-YQVRIQPSIHDHN Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade