

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09512 - IMMUNIZING PEPTIDE

SKU: EBP09512

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

 $NP_002947.1; \ NP_937883.1; \ NP_001234926.1; \ NP_001376220.1$ ID

 $\textbf{Immunogen} \begin{array}{l} \textbf{Peptide with sequence C-PHSTHHGSRGEERP, from the C Terminus of the protein sequence according to NP_002947.1; NP_937883.1; NP_001234926.1; NP_001376220.1. \end{array}$

The peptide was used in the production of EB09512 - Goat Anti-Restin / CLIP1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-PHSTHHGSRGEERP Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade