

EBP09513 - IMMUNIZING PEPTIDE

SKU: EBP09513



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_054706.1; NP_003364.1
Immunogen	Peptide with sequence C-RDYLDGSRGILSGT, from the internal regions (near N Terminus) of the protein sequence according to NP_054706.1; NP_003364.1.
Product	The peptide was used in the production of EB09513 - Goat Anti-Vinculin Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-RDYLDGSRGILSGT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png