

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09515 - IMMUNIZING PEPTIDE**

**SKU:** EBP09515

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

100 μg **Unit Size** 

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 001290.2 ID

Peptide with sequence C-KRQIFEPGLGMEH, from the internal region (near C Terminus) of the protein sequence Immunogen

according to NP 001290.2.

The peptide was used in the production of EB09515 - Goat Anti-Calponin 1 / CNN1 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-KRQIFEPGLGMEH Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade