

EBP09517 - IMMUNIZING PEPTIDE

SKU: EBP09517



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_008950.1; NP_861515.1; NP_861517.1
Immunogen	Peptide with sequence C-SGDKGISAFPESDN, from the internal region of the protein sequence according to NP_008950.1; NP_861515.1; NP_861517.1.
Product	The peptide was used in the production of EB09517 - Goat Anti-UBE2C / UBCH10 Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-SGDKGISAFPESDN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)