

EBP09534 - IMMUNIZING PEPTIDE

SKU: EBP09534



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006100.2; NP_001034708.1; NP_001269882.1; NP_001269883.1; NP_001269884.1; NP_001269885.1
Immunogen	Peptide with sequence RACREKDRDPEAQ, from the internal region of the protein sequence according to NP_006100.2; NP_001034708.1; NP_001269882.1; NP_001269883.1; NP_001269884.1; NP_001269885.1.
Product Comments	The peptide was used in the production of EB09534 - Goat Anti-PRMT5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	RACREKDRDPEAQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)