



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09534 - IMMUNIZING PEPTIDE

SKU: EBP09534

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 006100.2; NP 001034708.1; NP 001269882.1; NP 001269883.1; NP 001269884.1; NP 001269885.1 ID

Peptide with sequence RACREKDRDPEAQ, from the internal region of the protein sequence according to Immunogen

NP 006100.2; NP 001034708.1; NP 001269882.1; NP 001269883.1; NP 001269884.1; NP 001269885.1.

The peptide was used in the production of EB09534 - Goat Anti-PRMT5 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

RACREKDRDPEAQ Sequence

Shipping

Instructions Ambient

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png