

EBP09545 - IMMUNIZING PEPTIDE

SKU: EBP09545



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001129499.1; NP_006287.2; NP_001123954.1; NP_001123955.1
Immunogen	Peptide with sequence C-PNSQKVDSQKAATK, from the internal region of the protein sequence according to NP_001129499.1; NP_006287.2; NP_001123954.1; NP_001123955.1.
Product	The peptide was used in the production of EB09545 - Goat Anti-VRK2 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-PNSQKVDSQKAATK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)