

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09568 - IMMUNIZING PEPTIDE

SKU: EBP09568

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 004359.1; NP 958434.1 ID

Peptide with sequence C-DPGEVPEYPPYYQEE, from the internal region (near C Terminus) of the protein **Immunogen**

sequence according to NP 004359.1; NP 958434.1.

The peptide was used in the production of EB09568 - Goat Anti-Calponin 2 / CNN2 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-DPGEVPEYPPYYQEE Sequence

Shipping

Instructions Ambient

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade