

EBP09587 - IMMUNIZING PEPTIDE

SKU: EBP09587



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001838.2; NP_378667.1; NP_001274687.1; NP_001274688.1; NP_001274689.1
Immunogen	Peptide with sequence C-PSPEFETETLHNKES, from the internal region of the protein sequence according to NP_001838.2; NP_378667.1; NP_001274687.1; NP_001274688.1; NP_001274689.1.
Product Comments	The peptide was used in the production of EB09587 - Goat Anti-COL4A6 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PSPEFETETLHNKES
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)