

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09587 - IMMUNIZING PEPTIDE**

**SKU:** EBP09587

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 001838.2; NP 378667.1; NP 001274687.1; NP 001274688.1; NP 001274689.1 ID

Peptide with sequence C-PSPEFETETLHNKES, from the internal region of the protein sequence according to Immunogen

NP 001838.2; NP 378667.1; NP 001274687.1; NP 001274688.1; NP 001274689.1.

The peptide was used in the production of EB09587 - Goat Anti-COL4A6 Antibody. It is currently untested in **Product** 

Comments blocking applications. **Peptide** 

Sequence

C-PSPEFETETLHNKES

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png