

EBP09595 - IMMUNIZING PEPTIDE

SKU: EBP09595



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_002816.1
Immunogen	Peptide with sequence CGKLTGPKPQAESKK, from the internal region (near C Terminus) of the protein sequence according to NP_002816.1.
Product	The peptide was used in the production of EBP09595 - Goat Anti-Pleiotrophin / PTN Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	CGKLTGPKPQAESKK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png