

EBP09684 - IMMUNIZING PEPTIDE

SKU: EBP09684



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055423.1; NP_001311048.1; NP_001311053.1; NP_001311054.1; NP_001311055.1
Immunogen	Peptide with sequence C-TGSGRQEAQKTDAE, from the internal region of the protein sequence according to NP_055423.1; NP_001311048.1; NP_001311053.1; NP_001311054.1; NP_001311055.1.
Product Comments	The peptide was used in the production of EB09684 - Goat Anti-CYFIP1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TGSGRQEAQKTDAE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)