

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09684 - IMMUNIZING PEPTIDE

SKU: EBP09684

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

ID

Accession NP_055423.1; NP_001311048.1; NP_001311053.1; NP_001311054.1; NP_001311055.1

Peptide with sequence C-TGSGRQEAQKTDAE, from the internal region of the protein sequence according to

NP_055423.1; NP_001311048.1; NP_001311053.1; NP_001311054.1; NP_001311055.1.

The peptide was used in the production of EB09684 - Goat Anti-CYFIP1 Antibody. It is currently untested in Product

Comments blocking applications. Peptide

C-TGSGRQEAQKTDAE Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade