

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09704 - IMMUNIZING PEPTIDE

SKU: EBP09704

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_004407.2 ID

Peptide with sequence CHHIENVLKEDAR, from the internal region (near N Terminus) of the protein sequence

according to NP_004407.2.

The peptide was used in the production of EB09704 - Goat Anti-DTX1 / Deltex-1 Antibody. It is currently untested Product

Comments in blocking applications.

Peptide

CHHIENVLKEDAR Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade