

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09708 - IMMUNIZING PEPTIDE

SKU: EBP09708

SPECIFICATIONS

Formulation $100\mu g$ of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

 $NP_006051.1; NP_001278766.1; NP_001278767.1; NP_001207696.1; NP_001207697.1; NP_001278768.1; NP_001278768.1;$ Accession

NP_001207699.1; NP_001207700.1 ID

Peptide with sequence C-DTPYDSSASYEKEN, from the internal region of the protein sequence according to

Immunogen NP 006051.1; NP 001278766.1; NP 001278767.1; NP 001207696.1; NP 001207697.1; NP 001278768.1;

NP 001207699.1; NP 001207700.1.

The peptide was used in the production of EB09708 - Goat Anti-IKAROS / IKZF1 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-DTPYDSSASYEKEN Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade