

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09709 - IMMUNIZING PEPTIDE**

**SKU:** EBP09709

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_060243.2; NP\_001129662.1 ID

 $\label{eq:lower_lower_lower} \textbf{Immunogen} \ \ \frac{\text{Peptide with sequence NLTPSAHCINVRAS, from the internal region of the protein sequence according to NP_060243.2; NP_001129662.1.}$ 

The peptide was used in the production of EB09709 - Goat Anti-LAX1 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

**NLTPSAHCINVRAS** Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade