

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09744 - IMMUNIZING PEPTIDE**

**SKU:** EBP09744

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_000155.1; NP\_001295347.1 ID

Peptide with sequence C-EEQRQLPERAFR, from the internal region (near C-Terminus) of the protein sequence

according to NP\_000155.1; NP\_001295347.1.

The peptide was used in the production of EB09744 - Goat Anti-GIPR Antibody. It is currently untested in blocking Product

Comments applications.

Peptide

C-EEQRQLPERAFR Sequence

Ambient Instructions Shipping

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade