

EBP09752 - IMMUNIZING PEPTIDE

SKU: EBP09752



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	AAL47201.1; AAQ22368.1; EAX10059.1
Immunogen	Peptide with sequence CEHSSDFGWNGY, from the internal region of the protein sequence according to AAL47201.1; AAQ22368.1; EAX10059.1.
Product Comments	The peptide was used in the production of EBP09752 - Goat Anti-PRED6 / ABCCC13 Antibody. It is currently untested in blocking applications.
Peptide Sequence	CEHSSDFGWNGY
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)