

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09752 - IMMUNIZING PEPTIDE

SKU: EBP09752

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

AAL47201.1; AAQ22368.1; EAX10059.1 ID

Peptide with sequence CEHSSDFGWNGY, from the internal region of the protein sequence according to

AAL47201.1; AAQ22368.1; EAX10059.1.

The peptide was used in the production of EB09752 - Goat Anti-PRED6 / ABCCC13 Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

CEHSSDFGWNGY Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade