

EBP09807 - IMMUNIZING PEPTIDE

SKU: EBP09807



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001041666.1; NP_037414.3
Immunogen	Peptide with sequence NKRTAEQSCDQKLTN, from the internal region of the protein sequence according to NP_001041666.1; NP_037414.3.
Product Comments	The peptide was used in the production of EB09807 - Goat Anti-UHRF1 / ICBP90 Antibody. It is currently untested in blocking applications.
Peptide Sequence	NKRTAEQSCDQKLTN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)