

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09810 - IMMUNIZING PEPTIDE**

**SKU:** EBP09810

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

 $NP\_071910.3;\ NP\_001338019.1;\ NP\_001338020.1;\ NP\_001338021.1$ ID

Peptide with sequence C-RPTFIDRLANSLTKR, from the internal region of the protein sequence according to

Immunogen NP\_071910.3; NP\_001338019.1; NP\_001338020.1; NP\_001338021.1.

The peptide was used in the production of EB09810 - Goat Anti-IKZF4 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-RPTFIDRLANSLTKR Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade