

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP09821 - IMMUNIZING PEPTIDE

SKU: EBP09821

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_115807.1; NP_757385.1 ID

The peptide was used in the production of EB09821 - Goat Anti-SYVN1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-HERQHLEARLQS Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade