

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP09857 - IMMUNIZING PEPTIDE**

**SKU:** EBP09857

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Instructions Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_001089.1 ID

Immunogen Peptide with sequence C-QHVDVSPTSQRLQ, from the internal region of the protein sequence according to NP\_001089.1.

**Product** The peptide was used in the production of EB09857 - Goat Anti-Aconitase 2 Antibody. It is currently untested in

**Comments** blocking applications.

**Peptide** 

C-QHVDVSPTSQRLQ Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade