

EBP09879 - IMMUNIZING PEPTIDE

SKU: EBP09879



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Unit Size | 100 µg |
| Storage | Shipped at ambient temperature, store at -20°C. |
| Instructions | |
| Accession ID | NP_775298.1; NP_004662.2 |
| Immunogen | Peptide with sequence PQQVREICISRD, from the internal region of the protein sequence according to NP_775298.1; NP_004662.2. |
| Product Comments | The peptide was used in the production of EBP09879 - Goat Anti-PIAS2 (aa185-196) Antibody. It is currently untested in blocking applications. |
| Peptide Sequence | PQQVREICISRD |
| Shipping Instructions | Ambient |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |

DOCUMENTS

- [Data Sheet](#)