

EBP09884 - IMMUNIZING PEPTIDE

SKU: EBP09884



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004886.3; NP_899632.1; NP_001120933.1; NP_001120934.1
Immunogen	Peptide with sequence RDLYEKAKRDEPK-C, from the N-Terminus (near) of the protein sequence according to NP_004886.3; NP_899632.1; NP_001120933.1; NP_001120934.1.
Product Comments	The peptide was used in the production of EB09884 - Goat Anti-Cryopirin / NALP3 / NLRP3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	RDLYEKAKRDEPK-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)