

EBP09922 - IMMUNIZING PEPTIDE

SKU: EBP09922

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.**Unit Size** 100 µg**Storage Instructions** Shipped at ambient temperature, store at -20°C.**Accession ID** NP_076955.1**Immunogen** Peptide with sequence C-PTREYPSEQIIVTE, from the internal region of the protein sequence according to NP_076955.1.**Product Comments** The peptide was used in the production of EBP09922 - Goat Anti-DDA1 Antibody. It is currently untested in blocking applications.**Peptide Sequence** C-PTREYPSEQIIVTE**Shipping Instructions** Ambient**Product Grade** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)