



EBP09941 - IMMUNIZING PEPTIDE

SKU: EBP09941



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_848781.1
Immunogen	Peptide with sequence C-EERTINLPKSLKSH, from the C Terminus of the protein sequence according to NP_848781.1.
Product	The peptide was used in the production of EBP09941 - Goat Anti-Themis (mouse) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-EERTINLPKSLKSH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)