

EBP09997 - IMMUNIZING PEPTIDE

SKU: EBP09997



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_001727.1

Immunogen Peptide with sequence C-RESKSFTRSTVDT, from the C Terminus of the protein sequence according to NP_001727.1.

Product Comments The peptide was used in the production of EB09997 - Goat Anti-C5AR1/CD88 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-RESKSFTRSTVDT

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)