

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10027 - IMMUNIZING PEPTIDE

SKU: EBP10027

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_006227.1 ID

 $\label{eq:lmmunogen} \textbf{Immunogen} \begin{tabular}{l} \textbf{Peptide with sequence C-HMLSHAHQWVTAL}, from the internal region of the protein sequence according to NP_006227.1. \end{tabular}$

The peptide was used in the production of EB10027 - Goat Anti-POU3F3 / BRN1 / OCT8 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-HMLSHAHQWVTAL Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade