

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10048 - IMMUNIZING PEPTIDE

SKU: EBP10048

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

NP_001123645.1; NP_005366.2; NP_001123644.1; NP_001155128.1; NP_001155129.1; NP_001155130.1; Accession NP_001155131.1; NP_001155132.1 ID

Peptide with sequence C-HRDQTVKTQKENT, from the internal region of the protein sequence according to

Immunogen NP 001123645.1; NP 005366.2; NP 001123644.1; NP 001155128.1; NP 001155129.1; NP 001155130.1;

NP 001155131.1; NP 001155132.1.

The peptide was used in the production of EB10048 - Goat Anti-c-Myb (aa557-569) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-HRDQTVKTQKENT Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade