

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10056 - IMMUNIZING PEPTIDE

SKU: EBP10056

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_597677.2; NP_056199.2 ID

Immunogen Peptide with sequence C-EGNKSKVDTNKAHP, from the internal region of the protein sequence according to NP_597677.2; NP_056199.2.

The peptide was used in the production of EB10056 - Goat Anti-NIPBL (aa922-935) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-EGNKSKVDTNKAHP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade