

EBP10172-B - IMMUNIZING PEPTIDE

SKU: EBP10172-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001976.1; NP_001014763.1
Immunogen	Peptide with sequence C-PQRTAGVKVETTE., from the C Terminus of the protein sequence according to NP_001976.1; NP_001014763.1.
Product	The peptide was used in the production of EBP10172-B - Goat Anti-ETFB (aa 231-243), Biotinylated Antibody. It
Comments	is currently untested in blocking applications.
Peptide Sequence	C-PQRTAGVKVETTE.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)