

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10188 - IMMUNIZING PEPTIDE

SKU: EBP10188

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

 $NP_004997.4; \ NP_001186910.1; \ NP_001186911.1; \ NP_001186912.1; \ NP_001186913.1$ ID

 $\textbf{Immunogen} \begin{array}{l} \textbf{Peptide with sequence C-TEKSATYVNTEGR, from the internal region of the protein sequence according to NP_004997.4; NP_001186910.1; NP_001186911.1; NP_001186912.1; NP_001186913.1. \end{array}$

The peptide was used in the production of EB10188 - Goat Anti-NDUFS1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-TEKSATYVNTEGR Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade