

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10216 - IMMUNIZING PEPTIDE**

**SKU:** EBP10216

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_002337.1; NP\_032555.1 ID

Immunogen Peptide with sequence SCQEKDHYCIT, from the internal region of the protein sequence according to NP\_002337.1; NP\_032555.1.

**Product** The peptide was used in the production of EB10216 - Goat Anti-LY6E (mouse) Antibody. It is currently untested in

Comments blocking applications.

**Peptide** 

**SCQEKDHYCIT** Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade