

EBP10220 - IMMUNIZING PEPTIDE

SKU: EBP10220



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004172.2

Immunogen Peptide with sequence EVRFSAVALCK, from the C Terminus of the protein sequence according to NP_004172.2.

Product The peptide was used in the production of EB10220 - Goat Anti-UCHL1 (aa 211-221) Antibody. It is currently

Comments untested in blocking applications.

Peptide Sequence EVRFSAVALCK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)