

EBP10231 - IMMUNIZING PEPTIDE

SKU: EBP10231



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003816.2; NP_570710.1
Immunogen	Peptide with sequence C-TDKADTNRDR, from the C Terminus of the protein sequence according to NP_003816.2; NP_570710.1.
Product Comments	The peptide was used in the production of EBP10231 - Goat Anti-SNAP23 (aa 135-144) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TDKADTNRDR
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)