

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10235 - IMMUNIZING PEPTIDE

SKU: EBP10235

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_001271.1; NP_001287758.1; NP_001287744.1 ID

 $\label{eq:local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_local_$

The peptide was used in the production of EB10235 - Goat Anti-CIRBP (aa 81-91) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-QAGKSSDNRSR Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade