

EBP10258 - IMMUNIZING PEPTIDE

SKU: EBP10258



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005202.2; NP_001275634.1; NP_001362250.1
Immunogen	Peptide with sequence C-NDSFSEQDLDK, from the internal region of the protein sequence according to NP_005202.2; NP_001275634.1; NP_001362250.1.
Product	The peptide was used in the production of EBP10258 - Goat Anti-CSF1R (aa734-744) Antibody. It is currently
Comments	untested in blocking applications.
Peptide Sequence	C-NDSFSEQDLDK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)