

EBP10312 - IMMUNIZING PEPTIDE

SKU: EBP10312



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_083077.2

Immunogen Peptide with sequence CEKEDSSTCQ, from the internal region of the protein sequence according to NP_083077.2.

Product Comments The peptide was used in the production of EBP10312 - Goat Anti-Muc5b (mouse) Antibody. It is currently untested in blocking applications.

Peptide Sequence CEKEDSSTCQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)