

EBP10332 - IMMUNIZING PEPTIDE

SKU: EBP10332



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_004191.2

Immunogen Peptide with sequence C-DIPTEKLRET, from the internal region of the protein sequence according to NP_004191.2.

Product Comments The peptide was used in the production of EBP10332 - Goat Anti-Syt7 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-DIPTEKLRET

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)