

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10340 - IMMUNIZING PEPTIDE

SKU: EBP10340

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001139791.1; NP_076416.2; NP_001139792.1; NP_001139793.1; NP_001139794.1 ID

Peptide with sequence C-EAKPTRADAH, from the internal region of the protein sequence according to Immunogen

NP_001139791.1; NP_076416.2; NP_001139792.1; NP_001139793.1; NP_001139794.1.

The peptide was used in the production of EB10340 - Goat Anti-GRAMD3 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-EAKPTRADAH Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade