

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10358 - IMMUNIZING PEPTIDE

SKU: EBP10358

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

 $NP_001104753.1; \ NP_001104754.1; \ NP_001104755.1; \ NP_000609.1$ ID

Peptide with sequence C-RSVRAQRHTD, from the internal region of the protein sequence according to

Immunogen NP_001104753.1; NP_001104754.1; NP_001104755.1; NP_000609.1.

The peptide was used in the production of EB10358 - Goat Anti-IGF1 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-RSVRAQRHTD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade