

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP10371 - IMMUNIZING PEPTIDE**

**SKU:** EBP10371

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_057336.3 ID

Immunogen Peptide with sequence C-SREYDGNIRQAT, from the internal region of the protein sequence according to NP\_057336.3.

The peptide was used in the production of EB10371 - Goat Anti-BIRC6 / Apollon (aa 4724 -4735) Antibody. It is **Product** 

**Comments** currently untested in blocking applications.

**Peptide** 

C-SREYDGNIRQAT Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade