

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10383 - IMMUNIZING PEPTIDE

SKU: EBP10383

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_001010922.1 ID

 $\label{eq:lower_lower} \textbf{Immunogen} \ \ \ \ \ \ \ \ \ Peptide with sequence ESLSKTWCAQ, from the internal region of the protein sequence according to NP_001010922.1.$

The peptide was used in the production of EB10383 - Goat Anti-BCL2L15 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

Sequence

ESLSKTWCAQ

Shipping

Ambient Instructions

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png