

EBP10435-B - IMMUNIZING PEPTIDE

SKU: EBP10435-B



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_733936.1; NP_001071.1
Immunogen	Peptide with sequence C-SQDPAQIWRVAE., from the internal region of the protein sequence according to NP_733936.1; NP_001071.1.
Product	The peptide was used in the production of EBP10435-B - Goat Anti-ALDH5A1 (aa 485-496), Biotinylated Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-SQDPAQIWRVAE.
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)