

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP10446 - IMMUNIZING PEPTIDE

SKU: EBP10446

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_009165.1; NP_077324.1 ID

 $\label{eq:local_problem} \textbf{Immunogen} \ \ \ \ \ \ \ \ \ Peptide with sequence C-QVALGNISSKRE, from the internal region of the protein sequence according to NP_009165.1; NP_077324.1.$

The peptide was used in the production of EB10446 - Goat Anti-DCTN3 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

C-QVALGNISSKRE Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade