

EBP10558 - IMMUNIZING PEPTIDE

SKU: EBP10558



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_005268.1; NP_963885.1; NP_963886.1
Immunogen	Peptide with sequence C-QTYFEMARTAGD, from the internal region of the protein sequence according to NP_005268.1; NP_963885.1; NP_963886.1.
Product	The peptide was used in the production of EBP10558 - Goat Anti-GPM6A Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-QTYFEMARTAGD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)